**Para su mayor información, contacte a su representante local de Imeco**

*For further information, contact your Imeco local representative*

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*The excellency in quality and service is  
no longer a faraway horizon.*

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### The Result of the Evolution

After over six years of research and trials IMECO is the first company in the nation to release this new product to the market, introducing important improvements with respect to the products already available. In a joint venture with the INVAP, they have developed GOOD FIBER, ensuring a cementless acetabulum with higher durability, resistance and reliability, in response to the constant challenges presented, which require of cups that can guarantee:

- Higher resistance to the mechanical strain to which they are subjected.
- The correct insertion of the liner to the shell.
- An increase in the range of bone integration on the largest possible surface of the acetabulum.
- A greater proven biointegration.

GOOD FIBER, from IMECO, ensures biocompatibility of its components, with optimal mechanical qualities, and an excellent level of resistance to the various forms of corrosion.

Unlike other implants with similar characteristics, one of GOOD FIBER virtues is that it contains a polyethylene layer of ultra high molecular weight, thicker than those in existence, and its mechanism of blockage and fixation with pressure clips avoids movement, with which it provides a perfect match and integration between both components.

The characteristics of its fibrillar coating (Fiber Cover) guarantees the biological fixation on a larger surface of the acetabulum in the long run.

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### Fibrillar Coating - The evolution of fixation

GOOD FIBER presents a fibrillar coating, Fiber Cover, formed by titanium fibers synthesized with the **diffusion bonding** method. Such union is obtained through atom migration from one wire to another, process that **highly enhances bone integration**, guaranteeing both the initial and the definitive fixation for longer time. Previous studies on the coating, have resulted in a compact mesh, with porous cavities and high mechanical resistance.

The materials biocompatibility is guaranteed by the previous use of them both in general implants, as well as in hip prostheses.

The quality of the titanium used for its development - Ti 6 Al 4 V alloy 1, grade 5 titanium - guarantees a **higher resistance to corrosion and a lower level of toxicity**.


### Validation

With regard to its structure, some traits of similar types of pieces have been kept, which ensures the validity of the design of its components.

The shell manufacture is carried out under strict standards of control and quality, with the objective of meeting the dimensional demands of the matching surface, diminishing the tension of the contact and reducing polyethylene wear.

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With a minimum 7 mm, we achieve an increase in the resistance level of the liner, key piece in the useful life of the implant.

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The fibrillar coating Fiber Cover enhances bone integration, given its 1.5 mm thickness, the 51.8% of its surface covered by wire, and the remaining by interfibrillar pores of an average size of 280 microns.

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
### Structure. Load control, increasing the useful life

*In designing and developing the GOOD FIBER components, we took into consideration the close matching relationship between the acetabular component and the polyethylene liner.*

*Therefore, the processing of the shell is taken under strict tolerance, according to internal <sup>3</sup> processes and quality standards, in order to correctly place the polyethylene liner, respecting the congruency between both.*

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
### Polyethylene piece, thickness above the standard

*Taking into account that the polyethylene cup is under constant friction and tension with the femoral head, GOOD FIBER introduces an improvement, based in an increase in the liner thickness, taking it to a minimum 7 mm, above that which is used in other acetabula.*

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*The reduction in the shell thickness allows for a liner thickness above 7 mm for a nominal diameter of 52 and higher, using femoral heads of 28 mm. For nominal diameters of 48 and 50, 22 mm femoral heads should be used, in order to achieve higher than 7 mm liner thickness.*

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*Unlikely femoral necks, the strain at the acetabulum is relatively low. With a thickness of 3.5mm, the shell has a mechanical resistance higher than that it undergoes.*

*The manufacturing process allows to obtain diameters that vary no more than 0.1 mm with respect to the nominal ones.*

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*The contact surface is approximately 160°, depending on the acetabulum diameter*

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### Primary cup fixation

*In the images, we can observe two kinds of acetabula, identical except for the amount of holes for the fixation screws. Series 200 presents 2 and series 900, presents 9.*

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*Series 200  
The holes are 72° apart, and each at a 45° angle from the pole.*

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*Series 900  
Four of the holes are placed so that their axis forms a 30° angle with the pole, and are 90° apart. The remaining 5 holes form a 60° angle with the pole, and are 72° apart each.*

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### Characteristics of the holes

*The screw support has been dimensioned so that it resists the load it may transmit.*

*Studies indicate that the holes are prepared to sustain the highest possible tension, and though they might suffer some deformity, the screw is not dislodged.*

*Studies carried out after the primary fixation of the holes show that there is no lack of cohesion in the shell fibers, which explains the resistance and consistency of the mesh.*

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### Mechanical properties of the screws

*According to the results of the previously mentioned studies, although some deformity can be seen on the liner, it is located at the contact area between screw and hole, while the tension on the rest of the cup surface is low.*

*It is seen that under the load imposed, which is 68% higher than the maximum resistance of the bone, the acetabular cup cannot dislodge itself from the screw holding it.*

*Conclusions: The polar screw head lodging design is satisfactory from the structural point of view for the tolerance range proposed for the lip of that lodging.*

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*As can be seen the root diameter is 3 mm, while the top is 6.5 mm.*

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The liner is fixed to the shell by a pressure clips system. The edges of three pairs of clips at the rim of the shell are introduced in the liner slot. The clips possess a radial interference with the base of the slot, which permits for a firm fixation.

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#### Liner insertion and fixation – Exclusive system that prevents axial dislocation and rotation around the cup axis.

Its triple clip anchoring system optimizes the liner fixation to the implant on six supporting points, therefore avoiding dislocation and rotation, and diminishing the wear on the matching surface, with which the piece has a longer useful life, and the correct functioning is guaranteed for much longer. Its anchoring method diminishes the wear and the release of polyethylene particles from the interface between acetabulum and liner. As can be seen in the graphic, the clips indent the liner, which reduces the sustained deformity.

Liner expulsion trials have resulted in values higher than 80 kgs (reaching, in some cases, 200 kgs), which is higher than some acetabula in the market.

## GOOD FIBER - SURGICAL TECHNIQUE

## PREOPERATORY PLANNING

## Preoperative measuring with Acetabulum Templates

Good Fiber surgical technique has been carefully designed so as to guarantee an adequate fixation, and to ensure a biomechanically correct procedure. To be biomechanically correct, an arthroplasty requires that the center of the acetabulum be placed in quasi anatomical position.

Preoperative planning with complete standard radiology, both in the AP as in the lateral projections, should be used in order to visualize the ideal size and correct placing of the prosthesis. This can be adequately selected given the variety of options provided by the templates. It is important to standardize the exposures, to be able to compare x-rays in the post operative.

## Patient positioning and preparation

The Good Fiber acetabular component may be used with any of the standard approaches for Total Hip Arthroplasty, being the adequate exposure of the whole ring and bottom of the acetabulum the objective of the chosen approach. Below, we present a typical posterolateral approach.

## Acetabulum reaming

The patient should be positioned on his side, with the affected hip upwards. Both the position and the stability of the patient must be correct and secured prior to the placing of the surgical fields. To correctly place the acetabular component, the pelvis must remain stable and the leg should maintain mobility.

Once adequate exposure is achieved and any obstruction by osteophytes has been removed, the acetabulum must be prepared with a reamer, easily placed inside it. The external diameters available for the reamers range between 44 and 62 mm, increasing by 2 mm.

Before starting the procedure, the orientation of the reamer should be verified. The lateral inclination recommended is of 45° and the anteversion is of 10° to 30° approximately. Remember that the lateral positioning of the patient allows for the pelvis to be in a slight flexion of 15° to 20°.

Use progressively increasing reaming diameters, and remove the remaining articular cartilage and/or part of the subcondral plaque, until you reach partially bleeding bone tissue (2/3 of the bleeding acetabulum). There should remain as much of the subcondral plaque as possible. Then proceed to remove only the avascular sclerosed bone tissue. Do not completely remove the subcondral plaque (of a horseshoe shape) up to the depth of the acetabular pit, and never ream medially once the external plaque is pierced. If an adequately bleeding tissue cannot be achieved by reaming, use a higher diameter ream for a bigger peripheral reaming, without advancing medially (Fig. 1).

The two thirds of the acetabular ring should be healthy, in order to provide contention. Bleeding bone tissue can be more easily obtained in the posterior and superior parts, towards the ileon. To evaluate the strength and thickness of the posterior and anterior columns, these should be frequently palpated during the reaming.

To achieve an optimal fit, by pressure, the reaming should be 2 mm less than the acetabular component established by the template. A 54 mm cup requires a 2 mm reaming, while a 56 or higher cup, requires 2 mm of reaming. Additional fixation by screws could be needed.

With the Good Fiber system, a complete test reduction can be achieved, before impacting the acetabular component, which allows for a later correct placing and orientation. Use the dummy acetabulum with the cup size to implant.

The dummy acetabulum must be screwed at the end of the dummy cup, equal in size to the last reamer used, checking that the acetabular component covers the whole of the internal surface of the patient's acetabulum. Then the next dummy acetabulum must be placed, which is 2 mm larger than the last reamer used.

## IMPACTING THE ACETABULAR COMPONENT

The metallic Good Fiber acetabular component must be adequately screwed in the tap. By standard techniques, impact the acetabular component onto the clean and dry acetabulum (Fig. 4).

An adequate inclination should be maintained throughout the impaction (45°) and 20° to 30° of anteversion. Unscrew the tap once the component is impacted. For additional iliac fixation, the surgeon might choose to insert one or more 6.5 mm diameter cancellous bone screws. There are screws both for the 200 series as for the 900 series, 5

different lengths of self-screw cancellous bone screws, ranging from 20 to 40mm, in 5 mm increases.

Once the site for the screws is chosen, drill only with the drill-bushing through the holes in the acetabular component, using a flexible 3.2 mm drill.

To calculate the correct length of the screw use a depth measurer, previously removing the drill from it.

Extreme care should be taken not to drill the internal wall of the acetabulum so as not to harm the neurovascular pelvic structures. Using a cardanic screwdriver with hexagonal head, place the self-screwing screws completely.

To allow for the placing of the dummy polyethylene liner, first unscrew the Good Fiber tap. Place the dummy acetabulum, in the position previously established by the preoperative templates.

Then proceed to secure the adequate polyethylene liner (plastic liner) (fig. 3). For instance, a dummy acetabular cup of 56 mm will fit a plastic liner of 56 mm. The dummy acetabular component may be placed in any rotation position with the Good Fiber dummy polyethylene liner (plastic liner).

In tests with posterior rim, the bone adjacent to the indentation corresponding to the highest portion of the rim (posterior rim) should be marked with a bovie or methylene blue. This helps in correctly positioning the definitive implant.

The acetabulum line guide allows placing at 45° of lateral inclination and, according to the patient's requirements, anteversion of 0°, 10°, 20° or 30°.

## Implantation of the polyethylene liner (plastic liner)

The polyethylene component must be placed within the metallic acetabular component. During the test reduction, and when working with plastic liners of posterior rim, the indentation should be placed in the polyethylene rim, so that it matches the mark made at the bone rim of the acetabulum. Using a correctly sized polyethylene tap (22, 28 mm.) impact the polyethylene liner inside the acetabular component. Remember that the highest portion of the rim matches this indentation. As the polyethylene liner is impacted it will fit and fix itself within the titanium acetabular component.

Verify that the inner diameter of the liner matches the diameter of the femoral neck ball that will be used, because once it is impacted, it cannot be reused, and remember to place and orientate the definitive polyethylene component to the dummy component.

Next table shows the size options available; the inventory of polyethylene liners gives a range of options in diameters of 22 and 28 mm.

## IMPLANTATION OF THE POLYETHYLENE LINER

If it was necessary to change the polyethylene liner from an already implanted Good Fiber cup, use the polyethylene extractor.

Because of the damage resulting from the extraction of the polyethylene liner, it should be replaced. With the polyethylene extractor it is possible to remove the liner from an implanted Good Fiber cup.

## Post operative management

According to the individual and physical abilities of each patient, the post operative cares should be given, following standardized guidelines.

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## Primary fixation

Good Fiber has a uniform fibrillar layout which, together with the porous surface, greatly facilitates a quick primary fixation, consequently providing implant stability. The results of microscopic studies carried on the mesh show a proportionate surface distribution between the fibers and pores, therefore meeting the necessary conditions for achieving a bone integration process superior to the standard.

## High-compatibility porous surface

The objectives the cup must meet are: firstly to generate a larger surface of bone adhesion, and secondly to provide firm fixation for the liner. Fixation is influenced by the geometry of the matching surface and the type of coating chosen for bone integration.

Because of its characteristics, the fibrillar coating, fiber cover enhances bone

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

integration. Its thickness is 1.5 mm. and 51.8% of its surface consists in wire, while the rest is formed by the interfibrillary pores, which have an average size of 280 microns (optimal value within the accepted range between 40 and 400 microns), and which permit the development and placing of bone tissue.

The composition of the fibrillar coating is given by the union of the fibers, which is achieved at a high degree of pressure and vacuum; the fibers present (according to microscopic analysis) nearly zero porosity in its synthesis, therefore achieving higher resistance, without any contamination.

Fiber Cover is a porous coating formed by slim titanium fibers of commercial purity Grade 12, which is considered non toxic for the organism, and provides higher resistance to corrosion, apart from a lower rate of ions release.

With respect to its adhesion to the acetabulum, the adhesion trials carried out <sup>2</sup> show that because of the vacuum diffusion process used, it can tolerate up to 588 kgs/cm<sup>2</sup> (which is much higher than the 282 kgs/cm<sup>2</sup> found in other acetabula).

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280 microns is the average optimal pore size, because smaller pores obstruct bone growth and larger ones produce a growth which is slow and lacking in consistency.

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The high friction capacity of the porous coating Fiber Cover provides with a superior means for biological fixation by bone growth.

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