

SCREW DIAMETER (mm)								imeco [®] implantes ortopédicos		
Thread Diameter	1.5	2.0	2.7	3.5	4.0	4.5	5.5	6.5		
Screw Type	Cortex				Cancellous	Cortex	Cortex	Cancellous		
Catalog Number	200.8	201.8	202.8	204.8	---	214.8	---	---	---	
Self-tapping	200.0	201.0	202.0	204.0	206.0	207.0	214.0	212.0	216.0	
Non-self-tapping									217.0	
									218.0	
DRILL BIT and TAP DIAMETER (mm)								imeco [®] implantes ortopédicos		
Drill Bit for Gliding Hole	1.5	2.0	2.7	3.5	none	4.5	5.5	4.5		
	Paso 0,6	Paso 0,8	Paso 1,0	Paso 1,75	Paso 1,75	Paso 1,75		Paso 2,75		
GLIDING HOLE										
THREADED HOLE										
Drill Bit for Threaded Hole	1.1	1.5	2.0	2.5	2.5	3.2	4.0	3.2		
Tap	1.5	2.0	2.7	3.5	4.0	4.5	5.5	6.5		
Catalog Number	311.15	311.19 311.21 (long)	311.26 311.28 (long)	311.32 311.33 (long)	311.34	311.46 311.48	311.55	311.66		
Drive Type	1.5 mm/2.0 mm Cruciform			2.5 mm Hex		3.5 mm Hex				
									imeco [®] implantes ortopédicos	